AccuVein[®]

AccuVein AV400 User Manual

English	1
Español	17
Français	33
Deutsch	49
Dansk	65
Italiano	81
Nederlands	97
Português	113
Suomi	129
Svenska	145
Türkçe	161
Русский	177
اللغة العربية	193
한국어	209
繁體中文	225
ภาษาไทย	241

A health professional's guide for use
and operation of the AccuVein AV400

 AccuVein®

This AV400 User Manual and additional information can also be downloaded from www.accuvein.com

AV400 User Manual
© AccuVein Inc. 2013
AccuVein Inc.
40 Goose Hill Road
Phone: (816) 997-9400
Cold Spring Harbor, NY 11724
United States of America
www.accuvein.com

AV400 is manufactured for:
AccuVein Inc.
40 Goose Hill Road
Cold Spring Harbor, NY 11724
United States of America

EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands
Phone: +31.70.345.8570
Fax: +31.70.346.7299

Wavelength 642 nm/785 or 830 nm
Peak Power 45 mW / 25 mW
Pulse Energy 485 nJ / 270 nJ
Pulse Length 10.7 µs
IEC/EN 60825-1 (2001/07)

Device fully complies with
European Directive 2002/364/
EC. Conformité Européenne

For United States of America
customers: IEC/EN 60825-1 Safety
of Laser Product: Complies with 21
CFR 1040.10 and 1040.11 except for
deviations pursuant to Laser Notice
No. 50, dated June 24, 2007.

AccuVein® is the registered trademark of AccuVein Inc. The content of this manual including all photographs and drawings is the proprietary information of AccuVein Inc.

This manual last revised: October 2013

Warnings and Cautions

Read all instructions, cautions, and warnings prior to use. This product should be operated only by qualified medical professionals. The AV400 should not be used as the sole vein location method, and it is not a substitute for sound medical judgment and the visual and tactile location of veins.

Before using the AV400 on a patient, qualified medical professionals must read and understand this AV400 User's Manual. Before first use, users should compare how the AV400 detects veins with visual detection and palpation techniques.

English

AV400

AccuVein AV400 User Manual

A health professional's guide for use and
operation of the AccuVein AV400

English

AV400

AccuVein AV400 User Manual

A health professional's guide for use and operation of the AccuVein AV400

Basic Use and Operation	6	Charging, Cleaning and Troubleshooting	13	Warranty and Liability	14	Product Specification	15
Setup		Charging the AV400		Warranty Summary		Specifications	
Intended Use		Cleaning and Disinfecting		Limitation of Liability		Product labels	
Product Description		Troubleshooting					
User controls and screens							
Operating the AV400							

	Caution
	Consult operating instructions
	Refer to instruction manual/booklet
	The AV400 emits laser radiation
	Type B Applied Part
	Conformité Européenne (European Conformity). This symbol means that the device fully complies with European Directive 93/42/EEC
	ETL Classified means that most system elements conform to UL Std 60601-1
	Product Serial Number
	Reference to Catalog Model Number
	Universal Product Number
	European Representative
	Manufacturer

	European Union Waste Electrical and Electronic Equipment Directive Logo. Return is allowed for proper disposal.
	Temperature limitation
	Humidity limitation
	Do not use if package is damaged
	Keep dry
	Fragile, handle with care
	Reusable packaging
	Recyclable packaging
IPX0	No water ingress protection
Rx Only	US federal law restricts sale of this device only to physicians or other qualified medical professionals. For external use only.
	Direct Current. Example: 3.7 V indication of battery voltage
	Date of Manufacture

Warnings and Cautions

Patient and User Safety—Warnings

WARNING: The AV400 enables location of certain superficial veins and is not a substitute for sound medical judgment based on the visual and tactile location and assessment of veins. The AV400 should be used only as a supplement to the judgment of a qualified professional.

WARNING: The AV400 should only be operated when its battery is sufficiently charged (indicated when the battery icon at the upper right of the LCD screen is green), or when the device is operated in a powered hands-free accessory to ensure AV400 is available for use.

WARNING: AV400 vein location is dependent on a variety of patient factors and may not display veins on patients with deep veins, skin conditions, hair, scarring or other highly contoured skin surface, and adipose (fatty) tissue.

WARNING: The AV400 displays only superficial veins, and does so only to limited depths dependent on a variety of patient factors. The AV400 does not indicate vein depth.

WARNING: The AV400 emits Visible and Invisible laser radiation. Do not stare into beam. Class 2 laser product.

WARNING: Do not shine vein display light in eyes.

WARNING: Blink reflexes can be inhibited or reduced by disease, drugs, or other medical conditions. For patients who might have inhibited or reduced blink reflexes, AccuVein recommends outfitting the patient with a protective eye shield. A protective eye shield capable of blocking 642 nm and 785/830 nm should be used in this case.

WARNING: Stop using the AV400 if the red light does not turn on when the blue button is pushed.

WARNING: Do not hold the AV400 while performing venipuncture or other medical procedures.

WARNING: For external use only.

WARNING: Keep the AV400 and its battery out of the reach of children.

WARNING: To view vein location accurately, you must position the AV400 at the proper height, and directly over the center of the vein being assessed.

Patient and User Safety—Contraindications

WARNING: The AV400 should not be used to locate veins in or near eyes.

WARNING: The AV400 is not intended to be used as a diagnostic device or for treatment of any kind.

Patient and User Safety—Cautions

CAUTION: United States Federal law restricts this device to sale by or on the order of a physician or other qualified medical professional.

CAUTION: Operation or use of the AV400 in a manner different than specified in this AV400 User's Manual may result in hazardous radiation exposure.

CAUTION: Do not make any connections to the USB port on the AV400 device.

Equipment Care—Warnings

WARNING: Do not immerse the AV400 or the AV400 charging cradle in liquid or wet the AV400 such that liquid spills off.

WARNING: Do not attempt to open, disassemble, or service the battery pack. Do not crush, puncture, or dispose of the battery in fire or water. Do not short external contacts. Do not expose to temperature above 60°C / 140°F.

WARNING: Do not modify in any way the interior or exterior components of the AV400.

WARNING: To reduce the risk of fire or shock hazard and annoying interference, use only the recommended accessories and do not expose this equipment to rain or heavy moisture.

WARNING: Use only AccuVein accessories and replacement parts with the AV400. The use of non AccuVein accessories may degrade safety.

Equipment Care—Cautions

CAUTION: Use only AccuVein approved battery charging accessories, AccuVein BA300 and BA400 batteries, and AccuVein brand additional accessories. If the AV400 will not be used for an extended time, remove and store battery in a safe, dust-free location. The battery should be replaced only with an approved AccuVein battery. Removing or replacing the battery should not be done in the patient environment.

CAUTION: Electromagnetic Interference (EMI) can affect the proper performance of the device. Normal operation can be restored by removing the source of the interference.

CAUTION: Do not attempt to sterilize the AV400 with heat or pressure sterilization methods.

CAUTION: The AV400 will not display veins if operated outside its temperature range.

CAUTION: The AV400 may not display veins if operated under bright light such as bright sunlight.

CAUTION: Dispose of used batteries according to manufacturer's instructions. A recyclable lithium ion battery powers the AV400. Please call AccuVein at (816) 997-9400 for information on how to dispose or recycle this battery.

CAUTION: Do not disassemble or modify the AV400 or any of its charging accessories.

CAUTION: Do not service. The AV400 contains no customer serviceable components. The AV400 and its accessories should be serviced only by an authorized AccuVein repair department.

Basic Use and Operation

Setup

The AV400 Vein Illumination System includes an AV400 vein illuminator with rechargeable battery, a charging cradle, power supply and power adaptors, and instructional documents.

AccuVein recommends fully charging the AV400 before first use. Connect the charging cradle to the power supply. The AV400 comes with power supply country adaptors to ensure outlet compatibility in your country. If necessary, secure a different power adaptor to the power supply and before plugging it into an outlet.

Remove the plastic protective cover from the LCD before use

 **WARNING:** Inspect power cords and supplies for damage. If any damage is found do not use until replacement is obtained.

 **CAUTION:** Inspect the power adaptor on the power supply to confirm it is properly secured.

AccuVein also offers optional hands-free stands that many practitioners find useful (HF410, HF470). For information about these and other optional accessories visit the AccuVein web site www.accuvein.com.

Intended Use

The AccuVein® AV400 is a portable, hand-held instrument that helps medical professionals to locate certain superficial veins. The AV400 is intended to be used as a supplement to appropriate medical training and experience. The AV400 should not be used as the sole method for locating veins, and should be used only by a qualified medical professional, who should do so either prior to palpation to help identify the location of a vein, or afterwards to confirm or refute the perceived location of a vein. When using the AV400, medical practitioners should always follow the appropriate medical protocols and practices as required by their medical facility, as well as exercising sound medical judgment.

When used properly, the AV400 enables users to locate certain superficial veins in connection with medical procedures, such as venipuncture.

The AV400 can be used whenever the determination of vein location is appropriate.

Product Description

The AV400 operates by using infrared light to detect veins beneath the skin, then projecting the position of the veins on the skin surface directly above the veins. Qualified medical personnel can observe the vasculature as displayed to assist them in finding a vein of the right size and position for venipuncture and other medical procedures requiring the location of superficial veins. No training is required to operate the AV400.

The AV400 only shows superficial vasculature. The maximum depth that veins are displayed varies by patient. In addition, some patients' veins or a portion of their veins might not be displayed well or at all. Causes for less than optimal or lack of vein display include, but are not limited to, vein depth, skin conditions (e.g., eczema, tattoos), hair, scarring or other highly contoured skin surface, and adipose (fatty) tissue.


When held directly overhead, the AV400 accurately locates the center of a vein. Increasing the displacement from directly overhead results in an offset in the displayed vein position. Width of displayed vein may differ from the actual width depending on patient to patient differences and vein width. The center line of the vein is accurate when the device is being used correctly and should always be used as the target when performing venipuncture or other medical procedures.

The AV400 requires no routine or preventative maintenance.


The AV400 is portable internally powered by battery and approved for continuous operation. The AV400 is considered a Type B applied part and is not protected by ingress of water.

The technology in the AV400 is covered by the following patents: US 8,073,531; US 7,904,138; US 8,150,500; US 7,983,738; US 8,244,333; US 8,255,040; US 8,295,904; US 8,328,368; US 8,380,291; US 8,391,960; US 8,478,386; US 8,463,364; US 8,489,178; EP 1981395B1; Mexico 0357588. Other US and international patents pending. Visit www.accuvein.com for the most current information.

User controls and screens



Basic Use and Operation





AV400 controls at a glance	
	Power AV400 ON or OFF
	Turn vein display light ON or OFF
	Invert the projection
	About Product
	Select time-out

Display indications	
	Do not shine in eyes
	Battery charge status
	Next selection
	Select item
	Device in powered hands-free accessory


Basic Use and Operation

Operating screen display	
	Projection off, ready for use
	Projection on—dark veins over illuminated background
	Projection on—illuminated veins over dark background
	Selected time-out
	About product

Non-operating screen display	
	Unit too cold. Bring the device above 4°C/39°F
	Unit too hot. Bring the device below 33°C/90°F
	Less than 2 minutes of battery charge remaining. Charge battery immediately.
	Unit faulted. Call AccuVein Support at (816) 997-9400.

Operating the AV400

Positioning the AV400

Hold the device from 4–12" (10 to 30 cm) over the surface of the skin. Scan the area of interest to view . Once a vein is selected, make sure the vein display light is centered directly above the vein's center line. Don't move or tilt the device to either side of the vein or this will offset the projected vein from its true location beneath the skin. You can then often enhance display quality by slightly adjusting the height and angle to the skin. In particular, moving the device closer or further from the skin can help bring additional veins into view, depending on the patient's vasculature, room lighting, and depth of the veins.

⚠ WARNING: To view vein location accurately, you must position the AV400 directly over the center of the vein being assessed.

⚠ CAUTION: The AV400 projection may not display veins if operated under bright light such as bright sunlight

**Power on/off button**

Press the on/off button to power on or off the AV400.

**Blue button**

The blue button is located on the left side of the AV400. Press the blue button to power on the AV400 and turn on the vein display light. Press the button again to turn the vein display light off. This leaves the device on so you can navigate the device settings. Press the blue side button again to turn the vein display light back on and begin locating veins again.

**Time-out setting:**

The time-out setting determines the amount of time until the vein display light turns off. This setting may help increase battery operation time of the AV400. To set the time, press the left button to go into the time-out menu. The time out can be set to 0 (no time-out), 2 minutes, 5 minutes, or 15 minutes.


Operating the AV400 (continued)

Inverse setting:

The inverse setting inverts the dark (shadow) and light (red) projection areas of the AV400 such that the veins can be depicted by either light or dark. Changing the inverse setting may improve vein visibility and eye comfort level. AccuVein recommends reversing the vein display to show light veins in cases where veins are very small. When operating the device outdoors or in bright surroundings, use the dark veins settings to enhance visibility.

 **WARNING:** Do not shine vein display light in eyes.

 **WARNING:** Blink reflexes can be inhibited or reduced by disease, drugs, or other medical conditions. For patients who might have inhibited or reduced bright light blink reflexes, AccuVein recommends outfitting the patient with a protective eye shield. A protective eye shield capable of blocking 642 nm and 785/830 nm should be used.

Assess the Vasculature

While the vein display light is on, shine the AV400 over the patient's skin. You can do this before palpation, scanning quickly over the skin to help narrow down possible locations, or after palpation has confirmed vein location and suitability.

You can then often see veins better by rotating the AV400 slightly on its axis, and moving the device closer to or further away from the skin.

Confirm the Vein


After assessing the patient's vasculature, confirm the site for your procedure by verifying the location and suitability of the vein using normal medical techniques and good medical judgment, such as vein visualization, palpation, and other medical techniques.

If performing venipuncture

If you are performing venipuncture, you may use an AccuVein hands-free accessory so both your hands remain free to perform the procedure. You can also hand the AV400 to an assistant during the procedure. The AV400 should be used only by qualified medical professionals and only in conjunction with appropriate visual and tactile vein assessment techniques.


Charging the AV400

Verify the charging cradle is properly connected (see the Set-up section). Place the AV400 in the charging cradle with the LCD facing outward to begin charging. While the AV400 is in its charging cradle, the vein display light does not operate and the device cannot be used. An animated battery charge icon will display when the AV400 is charging. The light on the charging cradle will only turn on when a device is detected. If the light turns on when there is no device placed in the charging cradle please call AccuVein support for assistance.

To ensure the battery remains charged and ready for use, store the AV400 in its charging cradle or a powered AccuVein accessory.

CAUTION: Use only the AccuVein® PS310 power supply with the CC300 cradle and other AccuVein approved accessories.

When the AV400 battery gets low the battery indicator will flash yellow or red and beep repeatedly. Before the battery is completely depleted, an alert requesting to charge the battery will be displayed on the screen. The AV400 will resume normal operation once it has been charged.

In normal operation, the AV400 battery needs to be replaced about every two years. Longer charge cycles or shorter operating periods are indications that the battery should be replaced. To replace the battery, open the battery compartment on the narrow end of the device's underside. Use a screwdriver or edge of a small coin to unscrew the cover. Be careful not to touch or scratch the lenses on the back of the device. Do NOT use a tool to remove the battery.

Cleaning and Disinfecting

⚠ WARNING: Users should inspect the AV400 and clean and disinfect the AV400 as required by their institution's policies to ensure that it is sufficiently clean before each use.

⚠ WARNING: Do not use AV400 if the vein illumination window is scratched or dirty. The AV400 should be returned to AccuVein for servicing if the vein light window is scratched.

For peak performance, the optical surfaces and lenses on the back of the AV400 should be kept clean. Use alcohol wipes or any soft lens wipe with several drops of 70% isopropyl alcohol. Dirt or scratches on the vein display window show up as dark shadows in the vein projection.

To clean the body of the AV400 and its charging cradle, you may use a cloth moistened with 70% isopropyl alcohol, a 10% dilution of chlorine bleach in distilled water, or soap and water. You can also use common household disinfectants and wipes.


Do not immerse the AV400, the AV400 charging cradle, or the AV400 hands-free stand in liquid or get the AV400 or its components wet such that liquid spills off. Do not attempt to sterilize the AV400 with heat or pressure sterilization methods.

Do not clean the AV400 while in a charging cradle or hands-free accessory.

Do not clean the AV400 when the battery cover is removed.

Unplug the charging cradle before cleaning the cradle.

Troubleshooting

If the AV400 detects that it cannot operate properly, it turns off the vein display light and displays an alert or fault screen (See Non-Operating Display under User Screens). When the AV400 is not operating properly turn the device off using the round Power on/off button.

A yellow header on the screen indicates an alert. Alerts appear if you operate the device in an environment that is too cold (less than 4°C/39°F) or too hot (more than 33°C/90°F). Should a temperature alert occur, simply bring the AV400 back within its operating temperature range.

Other alerts indicate when the battery needs to be charged or replaced.

The fault screen has a red header and is picture above. If the fault screen appears, you should stop using the AV400 immediately. Turn off the device and contact AccuVein support via telephone at (816) 997-9400 option 3.

Warranty Summary

AccuVein warrants the AV400, when purchased new from an authorized seller, is of sound materials and workmanship and conforms to the manufacturer's specifications for a period of one year from the date of purchase from AccuVein or its authorized distributors or agents. The warranty period for accessories is as follows: Charging Cradle – 1 year; Power Supply – 1 year; Battery – 90 days. AccuVein will repair or replace any product that does not satisfy these warranties within this time limit at no cost to the customer. You must obtain a Return of Material Authorization (RMA) number from AccuVein support. The repaired unit has a 90-day warranty from the date it is received by the customer. Contact AccuVein support at (816) 997-9400 option 3 for service or additional information.

These warranties do not cover any misuse or abuse. These warranties are also void if the instrument is repaired by other than AccuVein or an authorized agent. There are no other warranties, either express or implied, other than those set forth in AccuVein's written warranty.

This is only a summary of the AccuVein warranty coverage. Contact AccuVein or see www.accuvein.com/legal/warranty for full details of warranty coverage.

Limitation of Liability

In no event shall either AccuVein or the customer be liable to the other for any incidental, indirect, special or consequential damages (including without limitation lost profits or lost business opportunity) that the other party may incur by reason of its having entered into or relied on a purchase or lease agreement, or arising out of the performance or breach of such an agreement, even if the party against whom a claim is made was advised or knew of the possibility of such damages. The foregoing limitation shall apply regardless of the form of the claim in which such liability may be asserted, including breach of contract, tort (including negligence) or otherwise.

Product Specification

Guidance and Manufacturer's Declaration			
The AV400, CC300, and HF470 are intended for use in the electromagnetic environment specified below. The customer or the user of the AV400, CC300, or HF470 should assure that it is used in such an environment.			
IMMUNITY Test	IEC 60601 Test Level	Compliance Level	Electromagnetic environment - guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst EC 61000-4-4	± 2 kV for power supply lines ± 1 kV input/output lines	± 2 kV for power supply lines Not applicable	Mains power quality should be that of a typical commercial or hospital environment
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_n (>95 % dip in U_n) for 0.5 cycles 40 % U_n (60 % dip in U_n) for 5 cycles 70 % U_n (30 % dip in U_n) for 25 cycles <5 % U_n (>95 % dip in U_n) for 5 s	<5 % U_n (>95 % dip in U_n) for 0.5 cycles 40 % U_n (60 % dip in U_n) for 5 cycles 70 % U_n (30 % dip in U_n) for 25 cycles <5 % U_n (>95 % dip in U_n) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the AV400 requires continued operation during power mains interruptions, it is recommended that the AV400 be powered from an uninterruptible power supply or a battery
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: U_n is the AC mains voltage prior to application of the test level.			

Recommended separation distances between portable and mobile RF communications equipment and the AV400			
The AV400, CC300, and HF470 are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the AV400, CC300, or HF470 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the AV400 as recommended below, according to the maximum output power of the communications equipment.			
Rated Maximum Output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.			

Guidance and Manufacturer's Declaration - Electromagnetic Emissions		
The AV400, CC300, and HF470 are intended for use in the electromagnetic environment specified below. The customer or the user of the AV400, CC300, or HF470 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic Environment Guidance
RF Emissions CISPR 11	Group 1	The AV400 uses RF energy only for its internal function. Therefore, its RF emissions are very low and not likely to cause any interference in nearby electronics
RF Emissions CISPR 11	Class A	The AV400 is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes
Harmonic Emissions IEC 6100-3-2	Complies	
Voltage Fluctuations IEC 6100-3-3	Complies	

Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
The AV400, CC300, and HF470 are intended for use in the electromagnetic environment specified below. The customer or the user of the AV400, CC300, or HF470 should assure that it is used in such an environment.			
IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the AV400, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 Vrms	$d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the AV400 is used exceeds the applicable RF compliance level above the AV400 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the AV400. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.


Product Specification

Specifications

Weight	275 grams (9.7 oz.)
Size	5 x 6 x 20 cm (2" x 2.4" x 7.9")
Battery	AccuVein® BA400 user replaceable Li-Ion 3.6V. BA400 – 3,100 mAh
Continuous (vein light on) run time on full charge:	Typically 180 minutes
Number of standard viewing procedures that can be performed using a fully charged battery	90
Charge time (full)	Typically 3 hours
CC300 + PS310	5V 2.0A 100V-240V 50Hz-60Hz 0.4A

Operating	Temperature	4°C to 33°C (39°F to 90°F)
	Humidity	5% to 85% RH non-condensing
	Pressure	75kPa to 106kPa
Transport	Temperature	-20°C to 50°C (-4°F to 122°F)
	Humidity	5% to 85% RH non-condensing
Storage	Temperature	-20°C to 50°C (-4°F to 122°F)
	Humidity	5% to 85% RH non-condensing

Product labels

Manufactured in USA for **AccuVein Inc.** **Rx ONLY. FOR EXTERNAL USE.**
Use BA400 battery. Conforms to UL Std. 60601-1; Cert. to CAN/CSA Std. C22.2 No 60111-M90
3.7 V **3100mAh** ETL CLASSIFIED

REF AV400 Rev. 1.0
SN AV12401234

0004742A

0004744A

Information Label with Laser Product Notice

Wavelength 642 nm / 785 nm
Peak Power 25 mW / 25 mW
Pulse Energy 270 nJ / 270 nJ
Pulse Length 11 µs
IEC / EN 60825-1 (2001/07)

EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands
Phone: +31.70.345.8570
Fax: +31.70.346.7299

0004742A

Laser Safety Label with Laser Classification Notice

VISIBLE AND INVISIBLE LASER RADIATION.
DO NOT STARE INTO BEAM. CLASS 2 LASER PRODUCT.

Complies with 21 CFR 1040.10 and 1040.11 except for deviations pursuant to Laser Notice No. 50, dated June 24, 2007.

AccuVein Inc.
40 Goose Hill Rd.
Cold Spring Harbor
New York 11724 USA
816.997.9400

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Guidance and Manufacturer's Declaration			
The AV400, CC300, and HF470 are intended for use in the electromagnetic environment specified below. The customer or the user of the AV400, CC300, or HF470 should assure that it is used in such an environment.			
IMMUNITY Test	IEC 60601 Test Level	Compliance Level	Electromagnetic environment - guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst EC 61000-4-4	± 2 kV for power supply lines ± 1 kV input/output lines	± 2 kV for power supply lines Not applicable	Mains power quality should be that of a typical commercial or hospital environment
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_n (>95 % dip in U_n) for 0.5 cycles 40 % U_n (60 % dip in U_n) for 5 cycles 70 % U_n (30 % dip in U_n) for 25 cycles <5 % U_n (>95 % dip in U_n) for 5 s	<5 % U_n (>95 % dip in U_n) for 0.5 cycles 40 % U_n (60 % dip in U_n) for 5 cycles 70 % U_n (30 % dip in U_n) for 25 cycles <5 % U_n (>95 % dip in U_n) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the AV400 requires continued operation during power mains interruptions, it is recommended that the AV400 be powered from an uninterruptible power supply or a battery
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: U_n is the AC mains voltage prior to application of the test level.

Recommended separation distances between portable and mobile RF communications equipment and the AV400			
The AV400, CC300, and HF470 are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the AV400, CC300, or HF470 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the AV400 as recommended below, according to the maximum output power of the communications equipment.			
Rated Maximum Output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

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
Manufactured in USA for **AccuVein Inc.** **Rx ONLY. FOR EXTERNAL USE.**
Use BA400 battery. Conforms to UL Std. 60601-1; Cert. to CAN/CSA Std. C22.2 No 601.1-M90
3.7 V 3100mAh ETL CLASSIFIED

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VISIBLE AND INVISIBLE LASER RADIATION. DO NOT STARE INTO BEAM. CLASS 2 LASER PRODUCT.

Complex with 21 CFR 1040.10 and 1040.11 except for deviations pursuant to Laser Notice No. 50, dated June 24, 2007.

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Guidance and Manufacturer's Declaration			
The AV400, CC300, and HF470 are intended for use in the electromagnetic environment specified below. The customer or the user of the AV400, CC300, or HF470 should assure that it is used in such an environment.			
IMMUNITY Test	IEC 60601 Test Level	Compliance Level	Electromagnetic environment - guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst EC 61000-4-4	± 2 kV for power supply lines ± 1 kV input/output lines	± 2 kV for power supply lines Not applicable	Mains power quality should be that of a typical commercial or hospital environment
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_n , (>95 % dip in U_n) for 0.5 cycles 40 % U_n , (60 % dip in U_n) for 5 cycles 70 % U_n , (30 % dip in U_n) for 25 cycles <5 % U_n , (>95 % dip in U_n) for 5 s	<5 % U_n , (>95 % dip in U_n) for 0.5 cycles 40 % U_n , (60 % dip in U_n) for 5 cycles 70 % U_n , (30 % dip in U_n) for 25 cycles <5 % U_n , (>95 % dip in U_n) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the AV400 requires continued operation during power mains interruptions, it is recommended that the AV400 be powered from an uninterruptible power supply or a battery
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: U_n is the AC mains voltage prior to application of the test level.			

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
Productetiketten

Manufactured in USA
for **AccuVein Inc.** Use BA400 battery.

Rx ONLY. FOR EXTERNAL USE.

3.7 V 3100mAh

Conforms to UL Std. 60601-1; Cert. to
CAN/CSA Std. C22.2 No 601.1-M90
ETL CLASSIFIED

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VISIBLE AND
INVISIBLE
LASER
RADIATION.
DO NOT STARE INTO
BEAM. CLASS 2
LASER PRODUCT.

Complies with 21 CFR 1040.10 and
1040.11 except for deviations
pursuant to Laser Notice No. 50,
dated June 24, 2007.

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Guidance and Manufacturer's Declaration			
The AV400, CC300, and HF470 are intended for use in the electromagnetic environment specified below. The customer or the user of the AV400, CC300, or HF470 should assure that it is used in such an environment.			
IMMUNITY Test	IEC 60601 Test Level	Compliance Level	Electromagnetic environment - guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV input/output lines	± 2 kV for power supply lines Not applicable	Mains power quality should be that of a typical commercial or hospital environment
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_n (>95 % dip in U_n) for 0.5 cycles 40 % U_n (60 % dip in U_n) for 5 cycles 70 % U_n (30 % dip in U_n) for 25 cycles <5 % U_n (>95 % dip in U_n) for 5 s	<5 % U_n (>95 % dip in U_n) for 0.5 cycles 40 % U_n (60 % dip in U_n) for 5 cycles 70 % U_n (30 % dip in U_n) for 25 cycles <5 % U_n (>95 % dip in U_n) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the AV400 requires continued operation during power mains interruptions, it is recommended that the AV400 be powered from an uninterruptible power supply or a battery
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: U_n is the AC mains voltage prior to application of the test level.			

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NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.
Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the AV400 is used exceeds the applicable RF compliance level above the AV400 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the AV400. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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
Manufactured in USA for **Accuvein Inc.** **Rx ONLY. FOR EXTERNAL USE.**
Use BA400 battery. Conforms to UL Std. 60601-1; Cert. to CAN/CSA Std. C22.2 No 601.1-M90
3.7 V 3100mAh ETL CLASSIFIED

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Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
The AV400, CC300, and HF470 are intended for use in the electromagnetic environment specified below. The customer or the user of the AV400, CC300, or HF470 should assure that it is used in such an environment.			
IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the AV400, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 Vrms	$d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.			
Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the AV400 is used exceeds the applicable RF compliance level above the AV400 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the AV400. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

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Manufactured in USA for **AccuVein Inc.** **Rx ONLY. FOR EXTERNAL USE.**
 Use BA400 battery. Conforms to UL Std. 60601-1; Cert. to CAN/CSA Std. C22.2 No 601.1-M90
 ETL CLASSIFIED

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Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
The AV400, CC300, and HF470 are intended for use in the electromagnetic environment specified below. The customer or the user of the AV400, CC300, or HF470 should assure that it is used in such an environment.			
IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the AV400, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 Vrms	$d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:
			
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.			
Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the AV400 is used exceeds the applicable RF compliance level above the AV400 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the AV400. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

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Manufactured in USA for **AccuVein Inc.** Use BA400 battery. Conforms to UL Std. 60601-1; Cert. to CAN/CSA Std. C22.2 No 801.1-M90

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Guidance and Manufacturer's Declaration			
The AV400, CC300, and HF470 are intended for use in the electromagnetic environment specified below. The customer or the user of the AV400, CC300, or HF470 should assure that it is used in such an environment.			
IMMUNITY Test	IEC 60601 Test Level	Compliance Level	Electromagnetic environment - guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst EC 61000-4-4	± 2 kV for power supply lines ± 1 kV input/output lines	± 2 kV for power supply lines Not applicable	Mains power quality should be that of a typical commercial or hospital environment
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_n (>95% dip in U_n) for 0.5 cycles 40% U_n (60% dip in U_n) for 5 cycles 70% U_n (30% dip in U_n) for 25 cycles <5% U_n (>95% dip in U_n) for 5 s	<5% U_n (>95% dip in U_n) for 0.5 cycles 40% U_n (60% dip in U_n) for 5 cycles 70% U_n (30% dip in U_n) for 25 cycles <5% U_n (>95% dip in U_n) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the AV400 requires continued operation during power mains interruptions, it is recommended that the AV400 be powered from an uninterruptible power supply or a battery
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: U_n is the AC mains voltage prior to application of the test level.			

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
Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the AV400 is used exceeds the applicable RF compliance level above the AV400 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the AV400. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Manufactured in USA for **AccuVein Inc.** **Rx ONLY. FOR EXTERNAL USE.**
 Use BA400 battery. Conforms to UL Std. 60601-1; Cert. to CAN/CSA Std. C22.2 No 601.1-M90
3.7 V 3100mAh ETL CLASSIFIED

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